

Bronopol

Catalog No: #S0880


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Bronopol |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Others |
| Calculated MW | 199.99 |
| Formulation | C3H6NO4Br |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Immunology
 SMILES: OCC(Br)(CO)[N+](=O)[O-]
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: 130-133°C
 Solubility: DMSO: 10 mM
 Appearance: White Crystal
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Matczuk M, et al. Int J Cosmet Sci. 2012 Oct;34(5):451-7.

Note: This product is for in vitro research use only