

6-Aminocaproic acid

Catalog No: #S0881


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	6-Aminocaproic acid
Brief Description	Inhibitors
Purification	99.88%
Target Name	Others inhibitor
Calculated MW	131.18
Formulation	C6H13NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular SystemSMILES:NCCCCC(O)=OPathways:OthersReceptor:PlasminogenBoling pt:Melting pt:207-209C?Solubility:Ethanol: 1 mg/mL (7.62 mM); Water: 26 mg/mL (198.21 mM); DMSO: 6 mg/mL (45.74 mM)Appearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Mochalkin I, et al. Biochemistry. 1999 Feb 16;38(7):1990-8.

Note: This product is for in vitro research use only