

D-Camphor


Catalog No: #S0882

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D-Camphor
Brief Description	Inhibitors
Purification	98.00%
Target Name	TRP/TRPV Channel inhibitor
Calculated MW	152.24
Formulation	C10H16O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C[C@@]12CCC(C1(C)C)CC2=O
Pathways: Membrane transporter/Ion channel
Receptor: TRPV3
Boiling pt: 204°C
Melting pt: 179-181°C
Solubility: soluble in Water
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Imming P, et al. Nat Rev Drug Discov. 2006 Oct; 5(10):821-34.

Note: This product is for in vitro research use only