

## Homotaurine


Catalog No: #S0883

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Homotaurine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Beta Amyloid inhibitor; GABA Receptor agonist
Calculated MW	139.17
Formulation	C3H9NO3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
SMILES: NCCCC(O)(=O)=O  
Pathways: Membrane transporter/Ion channel; Neuroscience  
Receptor: Aβ; GABA  
Boiling pt: Melting pt: 293°C  
Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
Stock solution can be stored below -20 °C for several months.

## References

1. Fariello RG, et al. Neurology, 1982, 32 (3): 241C5.

Note: This product is for in vitro research use only