

Behenic alcohol

Catalog No: #S0884


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Behenic alcohol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Antifection inhibitor
Calculated MW	326.6
Formulation	C22H46O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CCCCCCCCCCCCCCCCCCCCOPathways:Microbiology&virologyReceptor:GP340; GP350Boiling pt:180CMelting pt:65-72CSolubility:Ethanol: 4 mg/mL (12.24 mM)Appearance:White Granular solidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Leung DT, Sacks SL. Expert Opin PharmacOthers. 2004 Dec;5(12):2567-71.

Note: This product is for in vitro research use only