

4-Biphenylacetic acid


Catalog No: #S0888

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Biphenylacetic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX inhibitor
Calculated MW	212.24
Formulation	C ₁₄ H ₁₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system SMILES: OC(=O)CC1=CC=C(C=C1)C1=CC=CC=C1 Pathways: Immunology/Inflammation;

Neuroscience Receptor: COX Boiling pt: 365°C Melting pt: 159-160°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kohler C, et al. Arzneimittelforschung. 1980;30(4A):702-7.

Note: This product is for in vitro research use only