

Saccharin

Catalog No: #S0889


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Saccharin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	183.18
Formulation	C7H5NO3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:O=C1NS(=O)(=O)C2=CC=CC=C12Pathways:OthersReceptor:OthersBoiling pt:Melting pt:226-229CSolubility:Water: 4,000 mg/LAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Just T, et al. *Appetite*. 2008 Nov;51(3):622-7.

Note: This product is for in vitro research use only