

Furaltadone hydrochloride

Catalog No: #S0896


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Furaltadone hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	360.75
Formulation	C13H16N4O6HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:C1COCCN1CC2CN(C(=O)O2)/N=C/C3=CC=C(O3)[N+](=O)[O-].CIPathways:OthersReceptor:OthersBoiling pt:461.4CMelting pt:>300CSolubility:Water: 72 mg/mL (199.58 mM); DMSO: 17 mg/mL (47.12 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zolla L, et al. Biochem Cell Biol, 2005, 83(2), 166-175.

Note: This product is for in vitro research use only