

## Dehydroacetic acid

Catalog No: #S0899


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dehydroacetic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	168.15
Formulation	C <sub>8</sub> H <sub>8</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CC(=O)C1C(=O)OC(C)=CC1=O  
Pathways: Others  
Receptor: Others  
Boiling pt: 111-113°C  
Melting pt: 111-113°C  
Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Ohtsuki T, et al. Food Chem. 2013 Nov 15;141(2):1322-7.

Note: This product is for in vitro research use only