

Cetylpyridinium chloride monohydrate

Catalog No: #S0900


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cetylpyridinium chloride monohydrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	358
Formulation	C ₂₁ H ₃₈ ClNH ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: O.[Cl-].CCCCCCCCCCCCCCCC[N+]1=CC=CC=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: 81-84°C Solubility: Water: 50 g/L (at 20°C) Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Garcia-Godoy F, et al. Am J Dent. 2014 Aug;27(4):185-90.

Note: This product is for in vitro research use only