

## Chlorobutanol

Catalog No: #S0901


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Chlorobutanol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	177.45
Formulation	C <sub>4</sub> H <sub>7</sub> ClO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
 SMILES: O.CC(C)(O)C(Cl)(Cl)Cl.CC(C)(O)C(Cl)(Cl)Cl  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: 173°C  
 Melting pt: 78°C  
 Solubility: Soluble in DMSO  
 Appearance: White Crystal  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Chen SL, et al. Thromb Haemost. 1990 Nov 30;64(3):473-7.

Note: This product is for in vitro research use only