

Stugeron

Catalog No: #S0903


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Stugeron
Brief Description	Inhibitors
Purification	99.00%
Target Name	Calcium Channel inhibitor; Histamine Receptor antagonist
Calculated MW	368.51
Formulation	C26H28N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
 SMILES: C1CN(CCN1C/C=C/C2=CC=CC=C2)C(C3=CC=CC=C3)C4=CC=CC=C4
 Pathways: GPCR/G Protein; Membrane transporter/Ion channel; Immunology/Inflammation; Neuroscience
 Receptor: Calcium Channel; H1 receptor
 Boling pt: 509.2
 Melting pt: 117-120
 Solubility: Soluble in DMSO
 Appearance: White Crystal
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Terland O, et al. Neuropharmacology. 1999 Jun;38(6):879-82.

Note: This product is for in vitro research use only