

## Chloquinan(Chlorquinaldol)

Catalog No: #S0908


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Chloquinan(Chlorquinaldol)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	228.07
Formulation	C10H7Cl2NO
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:ImmunologySMILES:CC1=NC2=C(C=C1)C(Cl)=CC(Cl)=C2OPathways:OthersReceptor:OthersBoiling pt:350.7CMelting pt:108-112CSolubility:Ethanol: 6 mg/mL (26.3 mM); DMSO: 46 mg/mL (201.69 mM)Appearance:Gray PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Corrhons I, et al. Pathol Biol (Paris), 1991, 39(2), 136-139.

Note: This product is for in vitro research use only