

Metacetamol


Catalog No: #S0914

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Metacetamol |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Others |
| Calculated MW | 151.17 |
| Formulation | C ₈ H ₉ NO ₂ |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Nervous system
SMILES: CC(=O)NC1=CC(O)=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: 145-148°C
Melting pt: 145-148°C
Solubility: DMSO: 10 mM
Appearance: White Needle Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Hart A, Orr DL. *Antonie Van Leeuwenhoek*. 1975;41(3):239-47.

Note: This product is for in vitro research use only