

## Sulfaquinoxaline sodium

Catalog No: #S0917


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sulfaquinoxaline sodium
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	322.32
Formulation	C <sub>14</sub> H <sub>11</sub> N <sub>4</sub> NaO <sub>2</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: [Na+].NC1=CC=C(C=C1)S(=O)(=O)[N-]C1=NC2=CC=CC=C2N=C1 Pathways: Others Receptor: Others Boiling pt: 551.1°C Melting pt: >296°C Solubility: DMSO: 10 mM Appearance: Off-white Color Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. All NLM Curated PubMed Citations.

Note: This product is for in vitro research use only