

Trimebutine

Catalog No: #S0918


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trimebutine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Opioid Receptor agonist
Calculated MW	387.47
Formulation	C22H29NO5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CCC(COC(=O)C1=CC(OC)=C(OC)C(OC)=C1)(N(C)C)C1=CC=CC=C1Pathways:Endocrinology/Hormones; GPCR/G Protein; NeuroscienceReceptor:Opioid ReceptorBoiling pt:458CMelting pt:79CSolubility:DMSO: 50 mg/mL (129.04 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Prescire Int. 2013 Oct;22(142):241-2.

Note: This product is for in vitro research use only