

## Benfluorex hydrochloride

Catalog No: #S0919


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Benfluorex hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	387.82
Formulation	C19H21ClF3NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
 SMILES: Cl.CC(CC1=CC(=CC=C1)C(F)(F)F)NCCOC(=O)C1=CC=CC=C1  
 Pathways: Others  
 Receptor: HNF4α  
 Boling pt: Melting pt: 158-159°C  
 Solubility: Soluble in DMSO  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Lee SH, et al. ACS Chem Biol. 2013 Aug 16;8(8):1730-6.

Note: This product is for in vitro research use only