

Delta-Gluconolactone


Catalog No: #S0933

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Delta-Gluconolactone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	178.14
Formulation	C6H10O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: C([C@@H]1[C@H]([C@@H]([C@H](C(=O)O1)O)O)O)O Pathways: Others Receptor: LPH Boiling pt: 446.4°C Melting pt: 160°C Solubility: DMSO: 10 mM Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Pocker Y, Green E. J Am Chem Soc. 1973 Jan 10;95(1):113-9.

Note: This product is for in vitro research use only