

Cyclamic acid

Catalog No: #S0941


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Cyclamic acid |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Others |
| Calculated MW | 179.23 |
| Formulation | C6H13NO3S |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:Metabolism systemSMILES:OS(=O)(=O)NC1CCCCC1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:180CSolubility:Ethanol: 36 mg/mL (200.84 mM); Water: 36 mg/mL (200.84 mM); DMSO: 36 mg/mL (200.84 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Leban I, et al. Acta Crystallogr B. 2007 Jun;63(Pt 3):418-25.

Note: This product is for in vitro research use only