


## Description

Product Name	Amstat
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	157.21
Formulation	C8H15NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular SystemSMILES:NCC1CCC(CC1)C(O)=OPathways:OthersReceptor:PlasminogenBoiling pt:Melting pt:>300CSolubility:Water: 31 mg/mL (197.18 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Dunn CJ, Goa KL. Drugs. 1999 Jun, 57(6):1005-32.

Note: This product is for in vitro research use only