

Oxyclozanide

Catalog No: #S0955


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Oxyclozanide |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Others |
| Calculated MW | 401.46 |
| Formulation | C13H6Cl5NO3 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:ImmunologySMILES:C1=C(C=C(C(=C1Cl)O)NC(=O)C2=C(C(=CC(=C2Cl)Cl)Cl)O)ClPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sanabria R., et al. Efficacy of oxyclozanide against adult *Paramphistomum leydeni* in naturally infected sheep. *Vet Parasitol.* 2014 Dec 15;206(3-4):277-81.

Note: This product is for in vitro research use only