

Abamectin

Catalog No: #S0956


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Abamectin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	873.09
Formulation	C48H72O14
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CCC(C)C1OC2(CC3CC(CC=C(C)C(OC4CC(OC)C(OC5CC(OC)C(O)C(C)O5)C(C)O4)C(C)C=CC=C4COCC5C(O)C(C)=CC(C(=O)O3)C45O)O2)C=CC1 Pathways: Others Receptor: Antiparasitic Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Mogg TD, et al. Aust Vet J. 1990 Nov;67(11):399-401.

Note: This product is for in vitro research use only