

Ethopabate

Catalog No: #S0959


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ethopabate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	237.25
Formulation	C ₁₂ H ₁₅ NO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CCOC1=C(C=CC(NC(C)=O)=C1)C(=O)OC
 Pathways: Others
 Receptor: Folate metabolism
 Boiling pt: 148-151°C
 Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. El-Kosasy AM, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2015 Nov 5;150:430-9.

Note: This product is for in vitro research use only