

Iopanoic acid

Catalog No: #S0960


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Iopanoic acid
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	570.93
Formulation	C ₁₁ H ₁₂ I ₃ NO ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: CCC(CC1=C(I)C(N)=C(I)C=C1)C(=O)O
 Pathways: Others
 Receptor: Others
 Boiling pt: 153°C
 Melting pt: Solubility: Soluble in Water
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several months.

References

1. Renko K, et al. Endocrinology. 2012 May;153(5):2506-13.

Note: This product is for in vitro research use only