

## Trometamol

Catalog No: #S0961


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Trometamol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	121.14
Formulation	C4H11NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Respiratory systemSMILES:NC(CO)COPathways:OthersReceptor:OthersBoiling pt:167-172CMelting pt:2-8CSolubility:Water: 24 mg/mL (198.11 mM); DMSO: 2 mg/mL (16.5 mM)Appearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Holmdahl MH, et al. Acta Anaesthesiol Scand, 2000, 44(5), 524-527.

Note: This product is for in vitro research use only