

Histamine


Catalog No: #S0965

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histamine
Brief Description	Inhibitors
Purification	92.43%
Target Name	5-HT Receptor agonist
Calculated MW	111.15
Formulation	C5H9N3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: NCCC1=CN=CN1
Pathways: GPCR/G Protein; Neuroscience
Receptor: HT
Boiling pt: Melting
pt: Solubility: Free soluble in Water, Ethanol
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several months.

References

1. Hardie RC. Nature. 1989 Jun 29;339(6227):704-6.

Note: This product is for in vitro research use only