

Teicoplanin

Catalog No: #S0967


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Teicoplanin
Brief Description	Inhibitors
Purification	900ug/mg
Target Name	Others
Calculated MW	1709.39
Formulation	C77H77N5O31Cl2R
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CCCCCCCC(=O)N[C@@H]1[C@@H]([C@@H]([C@H](O[C@H]1OC2=C3C=C4C=C2OC5=C(C=C(C=C5)[C@H]([C@H]6C(=O)N[C@@H](C7=CC(=CC(=C7)C8=C(C=CC(=C8)[C@H](C(=O)N6)NC(=O)[C@@H]4NC(=O)[C@@H]9C1=CC(=CC(=C1)O)OC1=C(C=CC(=C1)[C@H](C(=O)N[C@@H](CC1=CC(=C(O3)C=C1)Cl)C(=O)N9)N)O)O)[C@@H]1[C@@H]([C@H]([C@@H]([C@H](O1)CO)O)O)O)C(=O)O)[C@H]1[C@@H]([C@H]([C@@H]([C@H](O1)CO)O)O)NC(=O)C)Cl)CO)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Water: 100 mg/mL (58.5 mM) Appearance: Off-white Color Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Berczki I, et al. Bioorg Med Chem Lett. 2014 Aug 1;24(15):3251-4.

Note: This product is for in vitro research use only