

Pyridoxine


Catalog No: #S0973

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyridoxine
Brief Description	Inhibitors
Purification	99.18%
Target Name	PLK
Calculated MW	169.18
Formulation	C8H11NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CC1=NC=C(CO)C(CO)=C1O
Pathways: Cell Cycle/Checkpoint
Receptor: PLK
Boiling pt: Melting pt: 214-215°C
Solubility: Soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Adams JB, et al. J Altern Complement Med. 2006 Jan-Feb;12(1):59-63.

Note: This product is for in vitro research use only