

Propoxur


Catalog No: #S0976

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Propoxur
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR agonist
Calculated MW	209.24
Formulation	C ₁₁ H ₁₅ NO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CNC(=O)OC1=CC=CC=C1OC(C)C Pathways: Neuroscience Receptor: AChR Boling pt: Melting pt: 91°C Solubility: Soluble in most polar solvents Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Mehta KD, et al. Environ Toxicol. 2014 Jun; 29(6): 705-13.

Note: This product is for in vitro research use only