

## Medrysone

Catalog No: #S0977


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Medrysone
Brief Description	Inhibitors
Purification	98.00%
Target Name	GR; Others
Calculated MW	344.5
Formulation	C22H32O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CC(=O)C12CC[C@H]3[C@@H]([C@@H]1CC[C@@H]2O)C(=O)C Pathways: Metabolism; Others Receptor: Glucocorticoid Receptor; Corticosteroid-binding globulin Boling pt: Melting pt: Solubility: Water: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Southren AL, et al. Invest Ophthalmol Vis Sci. 1983 Feb;24(2):147-52.

Note: This product is for in vitro research use only