

# (+)-Emetine dihydrochloride hydrate

Catalog No: #S0982


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	(+)-Emetine dihydrochloride hydrate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	571.58
Formulation	C <sub>29</sub> H <sub>40</sub> N <sub>2</sub> O <sub>4</sub> ·2HCl·H <sub>2</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

### Research

Area: Immunology SMILES: CC[C@H]1CN2CCC3=CC(=C(C=C3[C@@H]2C[C@@H]1C[C@@H]4C5=CC(=C(C=C5CCN4)OC)OC)OC)O.Cl.Cl.Pat  
 hways: Others Receptor: Others Boiling pt: Melting pt: 240-250°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Nomura T, Kutchan TM. J Biol Chem. 2010 Mar 5;285(10):7722-38.

Note: This product is for in vitro research use only