

## Artamine

Catalog No: #S0983


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Artamine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	149.21
Formulation	C5H11NO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
SMILES: CC(C)([C@H](C(=O)O)N)S  
Pathways: Others  
Receptor: Copper  
Boiling pt: 251.8°C  
Melting pt: 210°C  
Solubility: Water: 30 mg/mL (201.05 mM)  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Brewer GJ. Expert Opin PharmacOther. 2006 Feb;7(3):317-24.

Note: This product is for in vitro research use only