

## DAPK1 Polyclonal Antibody

Catalog No: #40832

Package Size: #40832-1 50ul #40832-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	DAPK1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	DAPK1 Polyclonal Antibody detects endogenous levels of DAPK1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human DAPK1 around the non-phosphorylation site of S308.
Target Name	DAPK1
Other Names	DAPK1; DAPK; Death-associated protein kinase 1; DAP kinase 1
Accession No.	Swiss-Prot: P53355NCBI Gene ID: 1612
Uniprot	P53355
GeneID	1612;
SDS-PAGE MW	160kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

## Images



Western Blot analysis of MOUSE-BRAIN cells using DAPK1 Polyclonal Antibody

---

Note: This product is for in vitro research use only