

DL-Penicillamine

Catalog No: #S0983L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DL-Penicillamine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	149.21
Formulation	C5H11NO2S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CC(C)(S)C(N)C(O)=OPathways:OthersReceptor:OthersBoiling pt:251.8CMelting pt:208-212CSolubility:Free soluble in Water.Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Brewer GJ. Expert Opin PharmacOthers. 2006, 7(3):317-24

Note: This product is for in vitro research use only