

Oxacillin sodium

Catalog No: #S0985


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Oxacillin sodium
Brief Description	Inhibitors
Purification	98.66%
Target Name	Others
Calculated MW	423.43
Formulation	C ₁₉ H ₁₈ N ₃ NaO ₅ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1=C(C(=NO1)C2=CC=CC=C2)C(=O)N[C@H]3[C@@H]4N(C3=O)[C@H](C(S4)(C)C)C(=O)[O-].[Na+] Pathways: Others Receptor: Others Boiling pt: 686.8°C Melting pt: 188°C Solubility: Water: 88 mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Woods GL, Yam P. Antimicrob Agents Chem Others. 1988 Nov;32(11):1614-8.

Note: This product is for in vitro research use only