

Gimeracil

Catalog No: #S0987


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gimeracil
Brief Description	Inhibitors
Purification	98.00%
Target Name	Autophagy inhibitor
Calculated MW	145.54
Formulation	C ₅ H ₄ ClNO ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: OC1=CC(=O)NC=C1Cl Pathways: Autophagy Receptor: DHPS Boiling pt: 281.9°C Melting pt: 274°C Solubility: DMSO: 29 mg/mL (199.25 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Prescire Int. 2013 May;22(138):122.

Note: This product is for in vitro research use only