

Desipramine hydrochloride

Catalog No: #S0991


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Desipramine hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor antagonist; Adrenergic Receptor
Calculated MW	302.84
Formulation	C18H23ClN2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: Cl.CNCCCNC1C=CC=CC=C2CCC2=CC=CC=C12
Pathways: Neuroscience; GPCR/G Protein Receptor: 5-HT; Adrenergic Receptor
Boiling pt: 407.4°C
Melting pt: 214-216°C
Solubility: Soluble in DMSO
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zavosh A, et al. Brain Res Bull. 1999 Jul 1;49(4):291-5.

Note: This product is for in vitro research use only