

Carbenoxolone disodium

Catalog No: #S0992


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carbenoxolone disodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Gap Junction Protein
Calculated MW	614.72
Formulation	C34H48Na2O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology
 SMILES: C[C@@]12CC[C@@]3(C(=CC(=O)[C@H]4[C@]3(CCC5[C@@]4(CC[C@@H](C5(C)C)OC(=O)CCC(=O)[O-])C)C)[C@@H]1CC(CC2)(C)C(=O)[O-])C.[Na+].[Na+]

Pathways: Cytoskeletal Signaling Receptor: Gap Junction Protein
 Boiling pt: 687.4°C
 Melting pt: Solubility: Water: 100 mM
 Appearance: White Powder

Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Hausmann W, Trnoky AL. Br J Pharmacol Chem Other. 1966 Feb;26(2):412-20.

Note: This product is for in vitro research use only