

Ethinodiol diacetate

Catalog No: #S0996


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ethinodiol diacetate
Brief Description	Inhibitors
Purification	99.48%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	384.51
Formulation	C24H32O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CC(=O)O[C@H]1CC[C@@H]2[C@H]3CC[C@]4([C@H]([C@@H]3CCC2=C1)CC[C@]4(C#C)OC(=O)C)CPathways:Endocrinology/HormonesReceptor:ERBoiling pt:461.1CMelting pt:128-131CSolubility:Water: 77 mg/mL (200.25 mM); DMSO: 77 mg/mL (200.25 mM)Appearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Siddique YH, et al. J Environ Biol, 2007, 28(2), 279-282.

Note: This product is for in vitro research use only