

Temephos

Catalog No: #S0997


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Temephos
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR agonist
Calculated MW	466.47
Formulation	C16H20O6P2S3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: COP(=S)(OC)OC1=CC=C(SC2=CC=C(OP(=S)(OC)OC)C=C2)C=C1 Pathways: Neuroscience Receptor: AChR Boling
 pt:518.5 Melting pt:30-31°C Solubility: Soluble in DMSO Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lima JB, et al. Am J Trop Med Hyg. 2003 Mar;68(3):329-33.

Note: This product is for in vitro research use only