

Dicloxacillin sodium

Catalog No: #S1001


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dicloxacillin sodium
Brief Description	Inhibitors
Purification	97.77%
Target Name	Others
Calculated MW	510.32
Formulation	C ₁₉ H ₁₈ Cl ₂ N ₃ NaO ₆ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1=C(C(=NO1)C2=C(C=CC=C2Cl)Cl)C(=O)N[C@H]3[C@@H]4N(C3=O)[C@H](C(S4)(C)C)C(=O)[O-].[Na+] Pathways:

Others Receptor: Others Boiling pt: 692.4°C Melting pt: 222-225°C Solubility: Water: 40mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Villegas-Guzman P, et al. J Environ Sci Health A Tox Hazard Subst Environ Eng. 2015;50(1):40-8.

Note: This product is for in vitro research use only