

Etodolac

Catalog No: #S1002


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Etodolac
Brief Description	Inhibitors
Purification	99.92%
Target Name	COX inhibitor; Retinoid Receptor inhibitor
Calculated MW	287.35
Formulation	C17H21NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CCC1=CC=CC2=C1NC1=C2CCOC1(CC)CC(O)=O
Pathways: Immunology/Inflammation; Metabolism; Neuroscience
Receptor: COX; RXR-alpha
Boiling pt: 507.9°C
Melting pt: 145-148°C
Solubility: Ethanol: 58 mg/mL (201.84 mM); DMSO: 58 mg/mL (201.84 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Glaser K, et al. Eur J Pharmacol. 1995 Jul 25;281(1):107-11.

Note: This product is for in vitro research use only