

# Proprantheline bromide


Catalog No: #S1006

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Proprantheline bromide
Brief Description	Inhibitors
Purification	95.00%
Target Name	AChR antagonist
Calculated MW	448.4
Formulation	C <sub>23</sub> H <sub>30</sub> BrNO <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism

system SMILES: [Br-].CC(C)[N+](C)(CCOC(=O)C1C2=CC=CC=C2OC2=CC=CC=C12)C(C)C

Pathways: Neuroscience Receptor: mAChR  
 Boiling pt: Melting pt: 159-161°C Solubility: Soluble in Water Appearance: White Solid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Saitoh H, et al. J Pharmacobiodyn. 1986 Dec; 9(12):1008-14.

Note: This product is for in vitro research use only