

## Aquacrine

Catalog No: #S1009


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Aquacrine
Brief Description	Inhibitors
Purification	99.09%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	270.37
Formulation	C18H22O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CCC2=O)CCC4=C3C=CC(=C4)O

Pathways:Endocrinology/HormonesReceptor:ERBoiling pt:Melting pt:Solubility:DMSO: 54 mg/mL (199.72 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Margeat E, et al. J Mol Biol. 2003 Feb 7;326(1):77-92.

Note: This product is for in vitro research use only