

## Vitamin E Acetate

Catalog No: #S1016


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Vitamin E Acetate
Brief Description	Inhibitors
Purification	99.12%
Target Name	Vitamin
Calculated MW	472.76
Formulation	C31H52O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous

SMILES: CC1=C2C(=C(C(=C1)OC(=O)C)CC[C@@](O2)(C)CCC[C@H](C)CCC[C@H](C)CCCC(C)C  
 Melting pt: 224°C Solubility: DMSO: 10 mM Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Beijersbergen van Henegouwen GM, et al. J Photochem Photobiol B. 1995 Jul;29(1):45-51.

Note: This product is for in vitro research use only