

Fluoroethyl


Catalog No: #S1017

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fluoroethyl
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	182.06
Formulation	C4H4F6O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:FC(F)(F)COCC(F)(F)FPathways:OthersReceptor:OthersBoiling pt:63.9CMelting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Rose L, Watson A. Anaesthesia. 1967 Jul;22(3):425-34.

Note: This product is for in vitro research use only