

Cefoxitin sodium

Catalog No: #S1023


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cefoxitin sodium
Brief Description	Inhibitors
Purification	96.71%
Target Name	Others
Calculated MW	449.43
Formulation	C16H16N3NaO7S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CO[C@@H]1([C@@H]2N(C(=O)C(=C(CS2)COC(=O)N)C(=O)[O-])NC(=O)CC3=CC=CS3.[Na+])Pathways: Others Receptor: Others Boiling pt: 843.4°C Melting pt: >160°C Solubility: DMSO: 32mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Salgado HR, Tozo GC. J AOAC Int. 2007 Mar-Apr;90(2):452-5.

Note: This product is for in vitro research use only