

## DDX52 Polyclonal Antibody

Catalog No: #40836

Package Size: #40836-1 50ul #40836-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	DDX52 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	DDX52 Polyclonal Antibody detects endogenous levels of DDX52 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human DDX52.
Target Name	DDX52
Other Names	DDX52; ROK1; HUSSY-19; Probable ATP-dependent RNA helicase DDX52; ATP-dependent RNA helicase ROK1-like; DEAD box protein 52
Accession No.	Swiss-Prot: Q9Y2R4NCBI Gene ID: 11056
Uniprot	Q9Y2R4
GeneID	11056;
SDS-PAGE MW	67kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

## Images



Western Blot analysis of HT29 cells using DDX52 Polyclonal Antibody

---

Note: This product is for in vitro research use only