

Moxonidine

Catalog No: #S1028


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Moxonidine
Brief Description	Inhibitors
Purification	98%+
Target Name	Others agonist
Calculated MW	241.7
Formulation	C ₉ H ₁₂ ClN ₅ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
 SMILES: COC1=C(NC2=NCCN2)C(Cl)=NC(C)=N1
 Pathways: Others
 Receptor: I1-imidazoline receptor
 Boiling pt: Melting pt: Solubility: Ethanol: 2 mg/mL (8.27 mM); DMSO: 24 mg/mL (99.3 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Aceros H, et al. Eur J Pharmacol. 2014 Oct 5;740:168-82.

Note: This product is for in vitro research use only