

Hesperidin

Catalog No: #S1035


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hesperidin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	610.56
Formulation	C28H34O15
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)OC[C@@H]2[C@H]([C@@H]([C@H]([C@@H](O2)OC3=CC(=C4C(=O)C[C@H](OC4=C3)C5=CC(=C(C=C5)OC)O)O)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:250-255C?Solubility:DMSO: 122 mg/mL (199.81 mM)Appearance:Gray PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Inderjit, Dakshini KM. J Chem Ecol. 1991 Aug;17(8):1585-91.

Note: This product is for in vitro research use only